

# Laboratory Distilled Water Plant

**LDWP**



SORAJ INSTRUMENTS Offers Laboratory Distilled Water Plant For Various R&D Sector, Laboratories, Educational Institute & Pharmaceutical Industries to get fast Distilled Water. Distillation removes inorganic solids, organics with higher boiling point than water, bacteria and pyrogens from water. Distillers can be turned on after standstill and instantly produce pure water without cleaning or replacement of filters. The distillation process offers autonomy and produces pure water without need for expensive product specific consumables like filters or cartridges. These water plants are engineered from premium quality raw materials, that are sourced from the authentic vendors of the market. The provided plant is suitable for different industries in order to get pure water for varied applications. Our offered plant is available for the clients in different capacities and specifications. This Distilled Water Plant is quality tested with the aid of diligent analysts.

Sl. No.	Parameter	Variants
1	Capacity	• 4.5Ltrs./hr., 10Ltrs./hr.
2	Auxiliary Supply	• 230V AC, 50Hz.
3	Purity	• Sterile 1 $\mu$ S/cm (1 MOhm*cm) @ 25 °C
4	Wattage	• 3KW & 6KW
5	Construction Material	• SS304 / SS316
6	Inlet Water Quality	• Tap Water
7	Mounting	• Wall

## Ordering Information For Model Selection:-

<b>LDWP</b>	1	2	3									
	<b>Capacity</b>	<b>Sheet Material</b>	<b>Mounting</b>									
	<table border="1"> <tr> <td>"4.5"</td> <td>For 4.5Ltr. / hr.</td> </tr> <tr> <td>"10"</td> <td>For 10Ltr. / hr.</td> </tr> </table>	"4.5"	For 4.5Ltr. / hr.	"10"	For 10Ltr. / hr.	<table border="1"> <tr> <td>"16"</td> <td>For SS-316</td> </tr> <tr> <td>"04"</td> <td>For SS-304</td> </tr> </table>	"16"	For SS-316	"04"	For SS-304	<table border="1"> <tr> <td>"W"</td> <td>For Wall Mounting</td> </tr> </table>	"W"
"4.5"	For 4.5Ltr. / hr.											
"10"	For 10Ltr. / hr.											
"16"	For SS-316											
"04"	For SS-304											
"W"	For Wall Mounting											

## Example:-

<b>LDWP</b>	4.5	16	W					
	<b>Capacity</b>	<b>Sheet Material</b>	<b>Mounting</b>					
	<table border="1"> <tr> <td>"4.5"</td> <td>For 4.5Ltr. / hr.</td> </tr> </table>	"4.5"	For 4.5Ltr. / hr.	<table border="1"> <tr> <td>"16"</td> <td>For SS-316</td> </tr> </table>	"16"	For SS-316	<table border="1"> <tr> <td>"W"</td> <td>For Wall Mounting</td> </tr> </table>	"W"
"4.5"	For 4.5Ltr. / hr.							
"16"	For SS-316							
"W"	For Wall Mounting							

\* Due to continuous development & improvements in design, we reserve the right to change the specification without notice



# SORAJ INSTRUMENTS & CONTROL ENTERPRISE

( An ISO 9001:2015 & 14001:2004 Certified & NSIC Approved Company )

Address: Shed-B1, Durgapur Industrial Estate, ( The WBSIDCL ), Phase-I, J.P.Avenue, Durgapur, West Bardhaman-713212

Contact: 09830469233, 09233392111, 8945531600

Email: sorajcontrol@gmail.com, admin@soraj.co.in, accounts@soraj.co.in

web: www.soraj.co.in